

ELECTRONIC SINGLE BREAST PUMP

MODEL NO: HNK-SUTPUMP-01

MADE IN
TURKEY



Hnimo®



Nimo ® Electronic Single Breast Pump is designed to be ergonomic and easy to use.

It has two different suction modes, bionic rhythm and classical fixed rhythm, which imitate the baby's suction rhythms. The device starts operating in low-level mode, helping to initiate the mother's milk flow. With the 4-stage suction level of both modes, the mother can choose the most appropriate rhythm for milking.

Thanks to the LED indicators, it can be easily seen which suction mode and which suction power it is in.

Not only with electric power but also it can be used with 4 AA batteries. For this reason it can be also use on travel or places where there is no electricity.

May Your Baby Always Be Well Fed

Milking with
adapter
or battery

Bionic rhythm
mimicking
suction

Painless,
comfortable
use

4 different
suction power
levels,
8 different
rhythm
options

Funnel
design that
allows to
manual using

Easy to use
with simple
user panel

150 ml
bottle
capacity

Automatic
shutdown
after 30
minutes

Quiet and
vibration-free
operation

Box Content of Product



1 piece
Nimo[®]
Electronic
Single
Breast Pump

1 Piece
DC/AC
Adapter

1 Piece
Carrying Bag

1 Piece
Feding Bottle
and Pacifier

1 Spare
Suction
Valve

1 Piece
Warranty
Card



EN 14350

**MADE IN
TURKEY**

**2 YEAR
WARRANTY**

Build Your Milk Storage For Your Baby



When your breasts feel full of milk and feel uncomfortable but your baby is not yet ready for feeding, you can relax and collect your milk without tiring yourself by using the Nimo[®] Electronic Single Breast Pump.

Greatest Companion of Moms



Don't worry about your milk being low. After breastfeeding your baby, milk your milk to increase your milk.

Comfortable Use with Adjustable Vacuum



Set the desired mode and vacuum level from the user panel and enjoy using the Nimo[®] Electronic Single Breast Pump.




Soft Silicone Chest Pillow

Thanks to the specially designed breast pillow, it facilitates the milk flow with massage feature while milking..

Milking with Adapter or Battery

You can easily milk with batteries in areas where there is no electricity such as during travels





%0 BPA


Nimo® Breast Pump
is made of BPA-free
materials

Nimo® Electronic Single Breast Pump has an inclined funnel design that will facilitate to expressing while sitting in an upright position and it has a breast cushion that will give comfort to your breast. It has a silicone diaphragm that prevents the expressed milk from escaping back to the hose and pump.

The parts that the milk touches, which is funnel, breast cushion, silicone diaphragm, silicone valve and feeding bottle are can be washe in the dishwasher.

- 1) Control and Motor Unit
- 2) AC/DC Adapter
- 3) Feeding Bottle
- 4) Single Way Milk Valve
- 5) Silicon breast pillow

- 6) Silicone expressing diaphragm
- 7) Body of Pump
- 8) Cover of Funnel
- 9) Air Hose
- 10) Connection of air hose and motor unit



1st Level
Indicator LED

2nd Level
Indicator LED

3rd Level
Indicator LED





4th Level
Indicator LED

Natural Rhythm
Suction Mode

Stable Rhythm
Suction Mode



Button Specifications

| Button | Function |
|---|-----------------------------|
|  | Switching device on and off |
|  | Change the milking mode |
|  | Increase the vacuum level |
|  | Decrease the vacuum level |

TECHNICAL SPECIFICATIONS

| | |
|-------------------------|---|
| Maximum Vacuum | 160 mmHg |
| Power | 6W-DC 6V 1A or (4x1.5 AA Alkaline Battery) (Batteries not included!) |
| Adapter | Input: 100-240V AC ~50/60Hz Output: DC 6V 1A |
| Automatic Shutdown Time | 30 minutes |
| Dimension | 105mm(Size) x 100mm(wide) x 60mm(Height) |
| Weight | 225 gr (Without Battery) |



New Generation
Technology