

# COMPRESSOR NEBULIZER

MODEL NO: HNK-NBL-S

MADE IN  
TURKEY



Hunkar®  
nimo

[www.hunkarmedikal.com](http://www.hunkarmedikal.com)



NIMO® Compressor Nebulizer is designed for use in the treatment of asthma, allergies and other respiratory ailments. NIMO® Compressor Nebulizer transforms the liquid medicine in the medicine container into aerosol vapor, making it easily inhalable thanks to the compressor motor that converts electrical energy into mechanical energy. NIMO® Compressor Nebulizer can give the medicine to the patient quickly and in the most appropriate way for the required treatment, without disturbing the patient. This device is suitable for use by patients of all ages. Thanks to its easy use and ergonomic structure, you can easily use it at home or at work.

Breathing  
is no longer  
hard thing  
thanks to  
**nimo<sup>®</sup>**

Lightweight and  
portable mini  
compact design

Low drug  
residue

Medicine  
container  
with adjustable  
steam outlet

High-density  
steam output  
thanks to the  
powerful  
compressor motor

Freedom of  
movement  
thanks to the  
long air hose

**FREE OF  
PHTHALATE  
&  
LATEX**

**2 YEAR  
WARRANTY**

## Medicine Container With Adjustable Steam Outlet

You can move the nebulization control valve located on the medication cup cover in the “+” or “-” direction. The amount of steam output is now under your control.



### + DIRECTION (Lower Airways)

- \*Cold
- \*Rhinitis
- \*Allergy Flu

### - DIRECTION (Upper Airways)

- \*Bronchitis
- \*BronchiolitisQ
- \*Asthma etc.

## Fast and Effective Treatment Opportunity

The size of the particles is between  $0.5 \mu\text{m}$  and  $5 \mu\text{m}$  and reaches even the lungs.

More intense vapours thanks to the powerful compressor motor

Comfort is with you thanks to NIMO®

Freedom of movement thanks to the long air hose

NIMO® is always right besides to you for comfortably breathes



NIMO<sup>®</sup> Compressor Nebulizer has been designed with your comfort with taking in our mind. It has an easy to use and ergonomic design.

We are always with you at your home or at your work thanks to its adjustable medicine container, long air hose, powerful motor etc.

## TECHNICAL SPECIFICATIONS


|                               |                 |
|-------------------------------|-----------------|
| Power                         | AC 220V 50Hz    |
| Power Consumption             | 150 VA          |
| MMAD                          | 3.3 $\mu$ m     |
| Air Outlet Range              | 4-6 L/dk        |
| Average Nebulization Rate     | 0.30 ml/dk      |
| Particle Size Range           | 0.5 - 9 $\mu$ m |
| Medicine Container's Capacity | 6 ml maks.      |
| Sound Level                   | < 60 dB         |



## BOX CONTENTS

|                             |                               |                  |                                   |
|-----------------------------|-------------------------------|------------------|-----------------------------------|
| 1 Nimo Compressor Nebulizer | 1 nozzle                      | 1 Adult Mask     | 1 Paediatric Mask                 |
| 1 Medical Container         | User Manual and Warranty Card | 1 Piece Air Hose | 5 Spare Air Filters for Nebulizer |



|  |   |   |
|--|---|---|
|  <p><b>T.C. SAĞLIK BAKANLIĞI</b><br/>Approved Product</p> |  <p><b>UTS</b><br/>ÜRÜN TAKİP SİSTEMİ</p> |  |
|   |    |  |

*All the best*  
**FOR YOUR HEALTH!**