

COMPRESSOR NEBULIZER

MODEL NO: HNK-NBL-S

MADE IN
TURKEY



Hunkar®
nimo

www.hunkarmedikal.com



NIMO® Compressor Nebulizer is designed for use in the treatment of asthma, allergies and other respiratory ailments. NIMO® Compressor Nebulizer transforms the liquid medicine in the medicine container into aerosol vapor, making it easily inhalable thanks to the compressor motor that converts electrical energy into mechanical energy. NIMO® Compressor Nebulizer can give the medicine to the patient quickly and in the most appropriate way for the required treatment, without disturbing the patient. This device is suitable for use by patients of all ages. Thanks to its easy use and ergonomic structure, you can easily use it at home or at work.

Breathing
is no longer
hard thing
thanks to
nimo[®]

Lightweight and
portable mini
compact design

Medicine
container
with adjustable
steam outlet

High-density
steam output
thanks to the
powerful
compressor motor

Freedom of
movement
thanks to the
long air hose

Low drug
residue

**FREE OF
PHTHALATE
&
LATEX**

**2 YEAR
WARRANTY**

Medicine Container With Adjustable Steam Outlet

You can move the nebulization control valve located on the medication cup cover in the “+” or “-” direction. The amount of steam output is now under your control.



+ DIRECTION (Lower Airways)

- *Cold
- *Rhinitis
- *Allergy Flu

- DIRECTION (Upper Airways)

- *Bronchitis
- *BronchiolitisQ
- *Asthma etc.

Fast and Effective Treatment Opportunity

The size of the particles is between $0.5\ \mu\text{m}$ and $5\ \mu\text{m}$ and reaches even the lungs.

More intense vapours thanks to the powerful compressor motor

Comfort is with you thanks to NIMO®

Freedom of movement thanks to the long air hose

NIMO® is always right besides to you for comfortably breathes



NIMO[®] Compressor Nebulizer has been designed with your comfort with taking in our mind. It has an easy to use and ergonomic design.

We are always with you at your home or at your work thanks to its adjustable medicine container, long air hose, powerful motor etc.

TECHNICAL SPECIFICATIONS

Power	AC 220V 50Hz
Power Consumption	150 VA
MMAD	3.3 μ m
Air Outlet Range	4-6 L/dk
Average Nebulization Rate	0.30 ml/dk
Particle Size Range	0.5 - 9 μ m
Medicine Container's Capacity	6 ml maks.
Sound Level	< 60 dB



BOX CONTENTS

1 Nimo Compressor Nebulizer	1 nozzle	1 Adult Mask	1 Paediatric Mask
1 Medical Container	User Manual and Warranty Card	1 Piece Air Hose	5 Spare Air Filters for Nebulizer



 <p>T.C. SAĞLIK BAKANLIĞI Approved Product</p>	 <p>ÜTS ÜRÜN TAKİP SİSTEMİ</p>	
		

All the best
FOR YOUR HEALTH!